RSC III® System

Reinforced Sternal Closure III® System

Application set for reinforcing the conventional closure of the sternum with wires following midline thoracotomy

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Sistema Reinforced Sternal Closure III

INSTRUCTION MANUAL

MANUEL D'INSTRUCTION

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RSC III® System

Reinforced Sternal Closure III® System

Description

The Reinforced Sternal Closure III® System is a method for reinforcing the conventional wire based closure of the sternum. A highly flexible, cone-shaped reinforcing element transfers the load from the sternal cerclage to a larger area of the sternum, thus preventing the wire from cutting into the cortical layers of the sternum. The complete lining of the bone channel in the sternum with the reinforcing element prevents also bleeding of the spongiosa to a large extent. This hemostatic effect is a further advantage to the conventional procedure.

Following midline thoracotomy for heart surgery or thoracic surgery where the sternum is bisected, a major post-operative complication is that of sternal instability. This can arise due to poor quality of the bone and post-operative respiratory movements such as excessive coughing.

The most important factor in preventing deep wound complications, chronic osteomyelitis, costal chondritis and mediastinitis is a stable sternal approximation. The Reinforced Sternal Closure III® System enhances significantly sternal stability.

The Reinforced Sternal Closure III® System is simple and rapid to use. It consists of wires without a needle, conic reinforcing elements and a single use thorn for the sternal punch. The sternal punch is a resterilizable surgical tool. The conic reinforcing elements, which are placed by the sternal punch directly through the sternum, transfer the load from the sternal cerclage to a larger area of the sternum. This prevents the wire from cutting into the cortical layers of the sternum. This first advantage increases the overall stability of the sternal osteosynthesis significantly. The second advantage is the additional hemostatic effect of the reinforcing element.

Indication

Patients at high risk of developing sternal instability following midline sternotomy.

Absolute Contraindication

Known allergy to stainless steel or its components

Relative Contraindication

Infection or suspected infection in the field of operation. In this case the surgeon must judge the possible advantages against the possible disadvantages or side effects.

Possible Side Effects

Possible specific side effects of the Reinforced Sternal Closure III® System may include:

- Acute and chronic foreign body reaction
- Allergy through hypersensitisation
- Chronic infection
- Injury to chest wall vessels
- Injury to the heart and/or bypasses
- Injury to the lung
- Migration of the reinforcing element

Contents and Packaging

The RSC III® application set is available in different specifications. Several wires with two reinforcing elements for each wire and one not implantable, single use thorn for the sternal punch are packaged in two sterile tubes and enclosed in an inside sterile peelable bag. Several bags are enclosed in a protective outer product carton. The tubes are closed on both ends with a strap and are transparent in order to identify its content.

The Reinforced Sternal Closure III® System fulfils Directive 93/42/EEC of the European Council for Medical Devices.

Sterilization

The RSC III® application set is sterilized by gamma irradiation. The contents of the bag remain sterile as long as the bag is undamaged, unopened and the directions concerning storage and expiry date are followed.

Warning: The RSC III® application set cannot be re-sterilized. If the above recommendations are ignored the user is solely liable for any undesirable consequences.

Storage

The RSC III® application set should be stored in the protective product carton in a cool dry place.

Expiry Date

The expiry date of the Reinforced Sternal Closure III® application set is indicated on the label of each product carton and envelope. It should not be used beyond the expiry date.

INSTRUCTIONS FOR THE USE OF THE RSC III® APPLICATION SET

General Remarks

The Reinforced Sternal Closure III® System may only be used by a surgeon. The RSC III® application set was manufactured by Alpha Research with great care and under stringent quality control. Since Alpha Research has no control over the indication and use of the device, Alpha Research cannot give further guarantees for the product. No representative of Alpha Research or distributor of the Reinforced Sternal Closure III® System is entitled to issue any additional guarantees.

Warning: General rules pertaining to asepsis must be observed when using the RSC III® application set. The reinforcing elements should only be used with the sternal closure wires supplied with each set. This prevents the possibility of electrocorrosion of either the reinforcing elements or the wires due to different material types.

Opening the Package

The bag should be peeled open and the sterile contents presented to the user under aseptic conditions. The sterile wires and reinforcing elements are then removed from the tube step by step by the user.

Implantation of the RSC III® System

First the single use thorn is removed from the tube and then screwed into the RSC® sternal punch. The instrument is opened and a reinforcing element slipped completely over the thorn. Then the sternal punch is placed rectangular to the midline sternotomy line and rectangular to the posterior surface of the sternum (Fig. 1).

Fig. 1



Fig. 2


The punch is closed and the reinforcing element pressed into the posterior cortical layer of the sternum. The thorn and the reinforcing element fully penetrate the sternum thus creating the channel for pushing through the corresponding wire at a later step (Fig. 2). The punch is then opened, removed from the body and loaded with the next reinforcing element.

Warning: In order to prevent loss in the operation site, the temporary placement of a compress under the sternum half before applying the procedure is recommended. Lost reinforcing elements must be removed as they can lead to uncontrolled damage of structures.

The punching and stapling procedure is repeated up to six or eight times at each half of the sternum. The posterior and anterior outlets of the bone channel are marked by the orifices of the reinforcing element. The sternal half is now lifted on one side. The wire is bent rectangular 3 cm from its end, grasped with a needle holder, pushed from posterior to anterior

through the reinforcing element orifices and the bone channel. The procedure is repeated at each placed reinforcing element on both corresponding sites thus preparing the single cerclages to be closed (Fig. 3).

Fig. 3

Warning: It is important to anchor the reinforcing elements completely into the bone. Partial anchoring at the edge of the sternum in the soft tissue will lead to a failure of the method.

After placement of the reinforcing elements and the wires, the wires are twisted together using the conventional technique. The risk of the wires cutting into the posterior and anterior cortical layers of the sternum is reduced by transferring the load from the sternum plates to a larger area of the sternum. Elevated tension can be achieved. This results in a significant increase of the stability of the sternal osteosynthesis. The complete lining of the bone channel with the reinforcing element prevents bleeding from the spongiosa to a large extent.

Warning: The RSC® sternal punch is to be closed with both hands under controlled application of pressure. The punching procedure is interrupted immediately after penetrating the ventral cortical layer in order not to crush the sternum und soft tissue. Penetration is to be noticed by a jerk of the RSC® sternal punch.

Disposal

The redundant reinforcing elements and wire offcuts should be put into a sharps bin for disposal.

RSC III® System

Reinforced Sternal Closure III® System

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Copyright© 2005
Alpha Research Deutschland GmbH